



**WIRELESS ULTRASOUND
PROBES
FOR MODERN
DIAGNOSTICS.**



CONTACT US

No.6/132 nagar nagar vayalur main road
Trichy-620017

Suloshanhealthcare@gmail.com
+91-9994248343, 9600199048

**SULOSHAN HEALTHCARE
& SERVICES**

The New Generation Medical
Tool



KEY HIGHLIGHTS



SPECIFICATIONS



- High-resolution imaging for breast lesion screening
 - Linear probe support (recommended for breast applications)
 - Color Doppler / PW Doppler capability
- Optional battery backup for mobility
 - 3-5 hours battery working time
 - Compatible with iOS, Android & Windows

ADVANCED WIRELESS ULTRASOUND TECHNOLOGY FOR MODERN HEALTHCARE

DISPLAY & SCANNING

Scanning mode: Electronic array.
Display: B, B/M, B+Color, B+PDI, B+PW. Probe elements: 80/128/192.
Frame rate: 18 fps. Image save: JPEG & AVI.

PROBE SPECIFICATIONS

Linear: 7.5/10MHz, 40mm head. Small Linear: 10/14MHz, 25mm head. Image Adjust: B-Gain, TGC, DYN, Focus, Depth, Harmonic, Color Gain.

ACCESSORIES

Bluetooth printer, wireless charger, needle puncture guide, tablet bracket, Android tablet/iPad, probe cover. Dimensions: 156×60×20mm, Weight: 220-250g.



M pro Interchangeable Probes – 10 Probe Options Available